

## Introduction

The functional network of human induced pluripotent stem cell (hiPSC)-derived neurons is a potentially powerful in vitro model for evaluating drug toxicity. Epileptiform activity is one of phenomena in neuronal toxicology. To evaluate the dynamics of epileptiform activities and the effect of anti-convulsant drug in cultured hiPSC-derived neurons, we used the high-throughput multielectrode array (MEA) system, where we simultaneously record extracellular potentials for 16 channels per well across 24-well plates. We firstly confirmed the modulation of activity by typical glutamatergic receptor antagonists/agonists in spontaneous firings. Spontaneous activities and pharmacological responses against synaptic related drugs were detected with high S/N ratio using high-throughput MEA system. Next, we examined chemically evoked epileptiform activity. Electrophysiological seizures were induced by pentylentetrazole (PTZ) and 4-Aminopyridine (4-AP), the most widely used chemical convulsant in animal models to screen for new anti-epilepsy drugs. We also examined the anti-convulsant effects of common clinical anti-epilepsy drugs (AEDs), phenytoin. PTZ and 4-AP induced a rapid increase in synchronized burst firings (SBFs) in a concentration-dependent manner. Phenytoin suppressed PTZ-induced epileptiform activity. From these results, we suggest that the electrophysiological assay in cultured human iPSC-derived neuron using high-throughput MEA system is a useful to investigate the neuronal toxicity in drug screening and pharmacological effects of human neurological disease.

## Material & Methods

### Human iPSC-derived neurons [XCell Science]

Human iPSC-derived neurons (XCell Science) were cultured  $3.0 \times 10^5$  cells/cm<sup>2</sup> on the MEA. Human iPSC-derived mature astrocyte (XCell Science) were added  $1.0 \times 10^5$  cells/well. After 8 days culture, medium was exchanged to BrainPhys medium (Stem cell technologies).

### High-Throughput MEA system [Alpha med scientific]

#### 24 wells (384 electrodes)


#### Recording



#### Low impedance and high sensitivity

#### Laster plots of spontaneous firings in hiPSC-derived neuros

Co-culture neuros with astrocytes (12 wells)      Neuros(12 wells)




## Conclusion


In conclusion, we detected epileptiform activities and effects of anti-epilepsy drugs in cultured hiPSC-derived neuronal networks and found that functional maturation at 6 weeks culture both neuros culture sample and co-culture neurons with astrocyte sample. High-throughput MEA system in cultured hiPSC-derived neurons proved useful for neuro pharmacological and neuro toxicological assays. Our results also provide an important indication for the toxicological evaluation using in vitro human neurons.

## Result 1

### Time course in spontaneous firings




## Result 2 Electrophysiological responses of NMDA and AMPA receptor at 6 weeks culture




## Result 3 Induction of epileptiform activity and effects of anti-epilepsy drugs

### 4-AP, Phenytoin



### PTZ, Phenytoin



The induction of epileptiform activity by 4-AP and PTZ and the suppressive effects by phenytoin were observed both neurons and co-culture with astrocyte samples.